




**First 2nd Generation
Swab Antigen Test**

**First Antigen Swab Test with a
sensitivity below 100 TCID₅₀/ml**

BioHit SARS-CoV-2 Antigen Rapid Tests (immuno-fluorescence) Nasopharyngeal Swab

Comparison of the Level of Detection of the most common antigen assays for SARS-CoV-2 from NP swab. Only the new test from BioHit deserves the attribute of a true 2nd Generation antigen test

Manufacturer	Name of Test	Limit of Detection** TCID ₅₀ /ml	Test Generation
BioHit*	SARS-CoV-2 Rapid Antigen	40	2nd
Siemens	CLINITEST® Rapid COVID-19 Ag Test	115	1st
Abbott	PanBio™ COVID-19 Ag Rapid	157	1st
Innova	SARS-CoV-2 Rapid Antigen	340	1st
Roche	SD Biosensor SARS-CoV-2 Antigen Rapid	494	1st
AlpStar	NADAL® COVID-19 Ag Rapid	502	1st

* The BioHit rapid assay also shows best in group results for clinical sensitivity (97,65%)

** Limit of detection is taken from the currently valid package inserts of each manufacturer

Visit us at www.dialane.org or write us at info@dialane.org