



Swab-based Anterior Nasal sampling

Simple, fast and painless Usage



BioHit SARS-CoV-2 Antigen Rapid Test for "Nasal Opening

- **Assured quality**
 - CE-IVD certified
 - BioHit is a TUV-South ISO13485 certified manufacturer
 - BfArM registered
- **Best in Class Performance**
 - 33 - 34 ct RT-PCR Equivalence
 - Independent evaluations
 - When it comes to COVID-19, demand the best
- **Analytical sensitivity. First 2nd generation test available on the market**
 - 50 pg/ml SARS-CoV-2 N antigen in saliva
 - Sensitivity: 96,1%
 - Specificity: 99,5%
- **Easy to use and safe for all**
 - Test kit contains all necessary materials (swabs, tubes, buffers, test cassettes)
 - **Very easy and quick, painless use**
 - Swab only 1 cm in the nostril and rotate approx. 8-10 times
 - Buffer contains detergents that inactivate viruses
 - Swab, swirl, shake, squeeze - result in 15 minutes
 - Disposal in household waste
- **Kit presentation**
 - 5-Kit / 25-Kit (1-Kit in work)